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PRINCIPLE OPERATION OF THE DEVICE:

The device "Omilon B-2" combines an automatic blood pressure meter and an indicator of blood glucose level. To determine the calculated concentration of glucose, it is necessary to measure the pressure on both arms. The method of measurement is based on the relationship between the dynamic elasticity of human vessels in the content of the hormone insulin and glucose in the blood. The device removes and analyzes the parameters of the pulse wave by a patented technique while measuring the pressure. Then the glucose concentration is mathematically calculated from these parameters.

DEVICE WORK PREPARATION:

It is necessary to install 4 batteries or a battery of size AA (finger) 1,2 - 1,5 V in the battery compartment, which is located at the bottom of the device. To do this, gently push the battery compartment cover and open according to the direction of the arrow. Be sure to observe the correct polarity while installing batteries! The polarity is indicated on the body of the device. A beep sound will confirm if the batteries are connected correctly. When the indication appears on the device, wait at least 2 minutes (until the indication disappears)! The device is ready for use. If the inflation time of the cuff exceeds 30 seconds, then the power supplies must be replaced.

The device switches off automatically in 2 minutes after the last press of any button! If the device is not used for a long time, it is necessary to remove the batteries from the device. Do not leave used batteries in the device. Push the buttons of the device without effort.

CUFF'S PREPARATION:

Pull the edge of the cuff approximately 5 cm into the metal ring. Put the cuff on your arm, while the tube should be pointed toward the palm of your hand. Wrap the cuff around the arm so that the lower edge of the cuff is 2-3 cm from the elbow fold. Secure the cuff so that it fits tightly around the arm, but does not squeeze it. Too tight or too loose cuffing can lead to inaccurate readings. If the arm is thick and has a pronounced taper, then it is recommended to put the cuff in a spiral. The device may show

inaccurate results if you wrap the sleeve of your clothes and squeeze your arm, preventing blood flow.

OPERATING PROCEDURE:

START - button to start / stop measurement;

SELECT - button to change the image information on the indicator or increase the pump pressure in the cuff during measurement;

MEMORY - the button for memorizing the measured parameters.

PRESSURE MEASUREMENT:

Do not eat, drink alcohol, smoke, exercise or take a bath before measuring the pressure. You must keep still while using an electronic device, so as to get correct reading before and during the measurement you need to calm down and relax, sit quietly for at least 3 minutes. The pressure is measured in the sitting position. Sit at the table, make yourself feel comfortable, put your feet on the floor and put your hand on the table. Press the "START" button to measure the pressure. The pressure is measured automatically: air is pumped into the cuff, then the device beeps at the same time when your heart beats, when the process is completed, air is released from the cuff. The display shows systolic ("upper") and diastolic ("lower") pressure at the end of the measurement.

Hold the "SELECT" button at the moment of inflation, if you have a pressure above 180 mm Hg! If it is necessary to stop the measurement, then press the "START" button. You must sit still quietly while measuring the pressure.

DETERMINATION OF THE CALCULATED CONCENTRATION OF GLUCOSE IN BLOOD:

The process of determining the calculated concentration of glucose in the blood with this device should be proceeded in the morning on an empty stomach or not earlier than 2 hours and 30 minutes after eating at any time of the day. To determine the calculated concentration of glucose, it is necessary to measure the pressure on both arms. First measure the pressure on the left arm. When the device is on, press the "START" button, the pressure is automatically measured. If the measurement is made correctly (there were no abrupt movements or conversation), it is necessary to press the "MEMORY" button, and the pressure results will be stored in the device memory. Then, measure the pressure on the right arm, but not more than 2 minutes after the measurement on the right arm, do not press the "MEMORY" button! After the measurement, you need to view the results by sequentially pressing and releasing the "SELECT" button, while setting the V mark according to the corresponding label on the front panel of the device.

The results are displayed in the following sequence:

- left arm pressure (systolic and diastolic pressure, V signs - no)*
- right arm pressure (systolic and diastolic pressure, V sign - below the right digit of diastolic pressure)*
- pulse (pulse result is in the right part of the indicator, sign V - at the bottom of the indicator, opposite the label - PULSE bpm)*
- the calculated concentration of glucose in mg / dl (sign V opposite the label - SUGAR / GLUCOSE mg / dl).*

- The calculated concentration of glucose in moles / liter (sign V opposite the label - SUGAR mmol / L).

REMEMBER: NEXT DEFINITIONS OF THE CALCULATED GLUCOSE CONCENTRATION OF THE SAME PATIENT DO NOT PROCEED LESS THAN 10 MINUTES AFTER THE LAST MEASUREMENT!

The device has two scales to display the glucose level:

1. Scale-1 (first) - is used for patients who do not have diabetes mellitus or patients with mild diabetes mellitus, whose glucose reduction to normal level is possible with diet compliance.

2. Scale-2 (second) - is used for patients with diabetes mellitus of medium severity who need sugar-lowering drugs (or who take hypoglycemic drugs) to normalize the level of sugar.

Press the buttons "SELECT" and "MEMORY" simultaneously if it's necessary to switch scales. When the second scale is on, the V mark appears in the lower right corner of the display.

The normal glucose concentration is 3.2 to 5.5 mmol / l or 60 to 100 mg / dl.