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If blanching is observed, i.e. the skin turns a white color, the injection should be stopped immediately and the area massaged until it returns to a normal color. After the injection, massage lightly the correction site to conform to the contour of the surrounding tissues. Retouching a correction undertaken with Cutegel should be immediately or a few weeks after the first injection.

WARNING

Cutegel is only intended for use as an intradermal implant.

Confirm that the product has not expired and that sterility has not been compromised before use, Product is for single use only; do not re-use.

In case of re-use, this can cause to depress Performance of the device and may lead to severe cross-infection

Used needles and syringes must be disposed into a collector designed for this purpose or in a biohazard sharps container

PRECAUTIONS FOR USE

- Cutegel is packaged for single patient use. Do not re-sterilize.
- Do not use if package is opened or damaged.
- Cutegel should not be injected into an area where there is an implant
- Cutegel should not be mixed with other products before implantation of the device.
- Hyaluronic acid products have a known incompatibility with quaternary ammonium salts such as benzalconium chloride
- Do not use in patients with bleeding disorders or patients who have undergone therapy with thrombolytic anticoagulants, or inhibitors of platelet aggregation in the preceding 2 weeks. In the same way, it is recommended to avoid taking aspirin, non-steroidal anti-inflammatory drug or high dose vitamin C the week before the injection
- Patients should be recommended not to apply make-up for 12 hours after the injection and to avoid prolonged exposure to sunlight, UV, as well as extreme cold and heat for two weeks after the injection.

ADVERSE EVENTS

Physicians must inform patients that there are potential adverse events that may occur or immediately be delayed after the injection.

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These adverse events include but are not limited to:

·Inflammatory reactions such as redness, swelling and tenderness may occur at the injection site. These reactions may last for two weeks.

·Nodules or induration are also possible at the injection site.

·Haematomas.

·Staining or discolouration of the injection site.

·Poor effect or weak filling effect.

·Cases of glabellar necrosis, abscess formation, granuloma and immediate or delayed hypersensitivity have been reported in the literature following hyaluronic acid injection. It is therefore important to take such possible complications into account.

Patients must report inflammatory reactions which persist for more than one week or any other secondary effect which develops to their practitioner as soon as possible. The practitioner should treat these as appropriate.

Any other undesirable side effects associated with injection of Cutegel must be reported to the distributor and/or to the manufacturer.

STERILIZATION INFORMATION SHELF LIFE AND STORAGE

Moist Heat Sterilization

Expiry date is indicated on each package.

2265 CE-marked according to MDD 93/42/EEC;
2265 is the Notified Body number for Cutegel

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